MARATHON™ XLPE Cemented Cup Hip Prosthesis

en-cs-et-hu-it-lv-pl-sk

Directory

en		Page 3
cs		Strana 12
et		Lehekülg 21
hu		31. oldal
lt		41 puslapis
lv		50 lpp.
pl		Strona 60
sk		Strana 71

MARATHON™ XLPE Cemented Cup Hip Prosthesis

IMPORTANT INFORMATION – Please Read Before Use

Some products may not be licensed in all jurisdictions.

STERILE

See product label for method

Instructions for use (IFU) website

For access to the e-IFU and any additional languages visit www.e-IFU.com

Implant card, patient information leaflet (PIL) and website

The implant card (and implant card label) is provided in the package with the implant. The information in the patient information leaflet (PIL) is based on the information in this IFU. The ePIL can be found on the website below.

For access to the ePIL visit www.ic.jnjmedicaldevices.com

DEVICE DESCRIPTION

This is an Ultra High Molecular Weight Polyethylene (UHMWPE) acetabular component with metal x-ray marker wires. A Total Hip Prosthesis is composed of individually packaged metal femoral hip stem, modular metal or ceramic femoral head and 2 piece metal-backed Ultra High Molecular Weight Polyethylene (UHMWPE), ceramic or all UHMWPE acetabular components designed to replace the natural articular surface of the hip joint.

MATERIALS AND RESTRICTED SUBSTANCES

The information below lists the substances in the implant. Additional trace substances may be present due to the manufacturing process. Please contact the address below for additional information.

These devices are a construct made from Ultra-High Molecular Weight Polyethylene (97-99% w/w) titanium alloy (1-2% w/w) and stainless steel alloy (less than 1% w/w) depending on the size.

Example of Acronym used on label(s):	Material Description:	Material nominally contains the following substances:
PE, UHMWPE, P1	Ultra-High Molecular Weight Polyethylene	Polyethylene (100%)
T1, Ti1	Titanium Alloy	Titanium (90%) Aluminium (6%) Vanadium (4%)
S2, SS2	Stainless Steel Alloy (316L)	Iron (62%) Chromium (18%) Nickel (14%) Molybdenum (3%) Manganese (2%)

HOW SUPPLIED

These device(s) are individually packaged and supplied STERILE. Remove from the package using accepted aseptic technique only after the correct size has been determined. Do not use these components if the sterile packaging appears damaged or unintentionally opened.

INTENDED USER

These devices are not to be used by unqualified personnel. It is essential that the surgeon and operating room staff are fully conversant with the appropriate surgical technique for the implant and associated instruments.

INTENDED USE (PATIENT TARGET GROUP)

Intended use

Total hip arthroplasty is intended to replace the damaged hip joint articulation in patients where there is evidence of sufficient sound bone to seat and support the components.

WARNING: Use only DePuy modular femoral heads and DePuy acetabular components with DePuy modular femoral stems. The outer diameter of the femoral head **MUST** be matched to the inner diameter of the acetabular component.

Cemented components

This acetabular device labelled "For cemented use only" is intended for use with bone cement.

PATIENT POPULATION

The intended patient population are the skeletally mature where there is evidence of sufficient sound bone to seat and support the components.

INDICATIONS

Total hip replacement is indicated in the following conditions:

1. A severely painful and/or arthritic joint (typically due to non-inflammatory degenerative joint disease)
2. Failed previous hip surgery
3. Acute fracture of the femoral head or neck

CONTRAINDICATIONS

The following conditions are contraindications for total hip replacement:

1. Active local or systemic infection.
2. Loss of musculature, neuromuscular compromise or vascular deficiency in the affected limb rendering the procedure unjustified.
3. Poor bone quality, such as osteoporosis, where, in the surgeon's opinion, there could be considerable migration of the prosthesis or a significant change of the fracture of the femoral shaft and/or the lack of adequate bone to support the implant(s).
4. Charcot's or Paget's disease.

WARNINGS AND PRECAUTIONS

CAUTION:

- Implants and trials components from different manufacturers or implant systems should never be used together.
- Hip prosthesis components should never be reimplanted. Even though the implant appears undamaged, the implant may have developed microscopic imperfections which could lead to failure.
- Always use a trial prosthesis for trial purposes. Trials should not be assembled with any components intended for permanent implantation. Trials must have the same configuration size, etc., as the corresponding components to be permanently implanted.
- Single use devices manufactured by DePuy have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Reuse can potentially compromise device performance and patient safety.
- Do not alter or modify implants in any way unless as described in the surgical technique.
- Quality bone stock is critical to ensure fixation of the device. Consider bone grafting and other reinforcement techniques, especially when performing revision surgery, to seat and support the device.

CAUTION: The following conditions, singularly or concurrently, tend to impose severe loading on the affected extremity, thereby placing the patient at higher risk for failure of the hip replacement:

1. Obesity or excessive patient weight.
2. Manual labour.
3. Moderate to high demand sports participation.
4. High levels of patient activity.
5. Likelihood of falls.
6. Alcohol or drug misuse addiction.
7. Other disabilities, as applicable.

CAUTION: The following conditions singularly or concurrently, tend to adversely affect the fixation of hip replacement implants:

1. Marked osteoporosis or poor bone stock.
2. Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g., diabetes mellitus, steroid therapies, immunosuppressive therapies, etc.).
3. History of general or local infections.
4. Severe deformities leading to impaired fixation or improper positioning of the implant.
5. Tumors of the supporting bone structures.
6. Allergic reactions to implant materials (e.g., bone cement, metal).
7. Tissue reactions to implant corrosion or implant wear debris.
8. Disabilities of other joints (i.e., knees and ankles).

NOTE: Diabetes, at present, has not been established as a contraindication. However, because of the increased risk for complications such as infection, slow wound healing, etc., the physician should carefully consider the advisability of hip replacement in the severely diabetic patient.

The surgical and postoperative management of the patient must be carried out with due consideration for all existing conditions. Mental attitudes or disorders resulting in a patient's failure to adhere to the surgeon's orders may delay postoperative recovery and/or increase the risk of adverse effects including implant or implant fixation failure.

Excessive physical activity or trauma to the replaced joint may contribute to premature failure of the hip replacement by causing a change in position, fracture, and/or wear of the implants.

EXPECTED CLINICAL BENEFITS

Total hip arthroplasty is expected to provide increased patient mobility and reduce pain by replacing the damaged hip joint articulation in skeletally mature patients where there is evidence of sufficient sound bone to seat and support the components.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

Some countries or regions have hip replacement registries which monitor the performance of hip implants. This registry data was analysed and shows that for a person's first total hip operation generally 95% (95 out of 100) last more than 10 years. This means that at 10 years about 5% (5 out of 100) of patients may have had additional operations to remove part(s). For revision hip replacement generally 80% last more than 10 years. This device meets this general implant lifetime. The risk of revision for an individual patient will depend on patient specific factors such as their anatomy, medical condition and lifestyle (e.g. weight and activity level), the surgical procedure and process.

A device specific SSCP can be found at the following link:

<http://ec.europa.eu/tools/eudamed>

Note: The Eudamed link will only be available after the Eudamed database goes live.

Reporting of Serious Incidents

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- a) the death of a patient, user or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat.

INFORMATION FOR USE

Preoperative

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY AND THE PROSTHESES WITH THE PATIENT BEFORE SURGERY. The discussion should include the limitations and possible consequences of joint replacement, the necessity to follow the surgeon's instructions postoperatively and the information provided in the patient leaflet.


The preoperative planning and surgical techniques for implantation of these hip replacement components evolved from the surgical experience gained during the development of many hip prostheses. Surgeons should not begin the clinical use of any hip prosthesis before they have thoroughly familiarized themselves with its specific implantation technique. Surgical technique brochures are available from DePuy.

Intraoperative


It is recommended that components at least one size larger and one size smaller than were preoperatively determined be available at surgery to accommodate intraoperative selection of the appropriate size(s).

Protective covers should be left on until the components are ready to be implanted. Do not use components if they have been dropped or have impacted a hard surface. Damage to the component may not be visible but could cause early failure of the prosthesis. The surgical instruments identified in the surgical technique have been designed for this arthroplasty procedure and must be used as shown in the surgical technique brochure.


When performing revision surgery, it is imperative that any fibrous membrane encapsulating the failed prosthesis be removed in its entirety. Residual fibrous membranes may interfere with implant fixation. Review any pre-existing devices for damage to mating surfaces that would prevent adequate seating or locking of the new device. Improper selection, placement, or positioning of implants may result in unusual stress conditions and a subsequent reduction in the functional life of the implant. Increasing the anteversion of an acetabular component in total hip replacement may result in instability and/or dislocation. Prior to closure, the surgical site should be thoroughly cleansed of bone chips, extraneous bone cement, ectopic bone, etc. Foreign particles at the metal/plastic or ceramic/plastic interface may cause excessive wear. Range of motion should be thoroughly checked for improper mating, instability, or impingement and corrected as appropriate.

1.

2.

3.

X-Ray Marker Wire Assembly

Postoperative

Strict adherence by the patient to the surgeon's instructions and warnings is extremely important. Accepted practices should be followed in postoperative care.

The patient should be released from the hospital with complete written instructions and warnings regarding exercises and therapies and any limitations on their activities. Component malposition or the use of constrained hip or total hip with modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

A continuing periodic patient follow-up is recommended because of the unknown functional lifetime of the implant in the specific patient, particularly with respect to the maintenance of implant fixation and bearing surfaces. Radiological investigations should be performed at each follow-up and compared with previous radiological investigations and correlated with the clinical assessment of the patient. If any radiological changes are observed, such as the occurrence of radiolucencies, bone resorption, or any changes in the position of an implant, these changes should be closely monitored to determine whether they are static or progressive and the patient treated appropriately.

ADVERSE EVENTS AND COMPLICATIONS

The following are generally the most frequently encountered adverse events and complications in hip arthroplasty:

General

The incidence and severity of complications in hip replacement are usually greater in revisions than in primary operations. Common problems encountered in revision surgery may include difficulty in placement of the incision, removal of ectopic bone and old bone cement, positioning and fixation of components, and/or obtaining adequate bony support. In general, increased operative time, blood loss, and risk of infection, pulmonary embolus and wound hematoma can be expected with revision procedures.

General

1. Change in position of the prosthetic components, often related to the factors listed in WARNINGS AND PRECAUTIONS.
2. Early or late loosening of the prosthetic components, often related to factors listed in WARNINGS AND PRECAUTIONS.

3. Fatigue fracture of the femoral stem or cup component, often related to factors listed in WARNINGS AND PRECAUTIONS.
4. Excessive wear or fracture of the bearing components due to: intraoperative damage to the prosthetic components, loose cement, bone fragments, metallic particles, ceramic particles or other factors listed in WARNINGS AND PRECAUTIONS.
5. Early or late infection.
6. Peripheral neuropathies. Subclinical nerve damage may also occur as a result of surgical trauma.
7. Tissue reactions, osteolysis, and/or implant loosening caused by metallic corrosion, allergic reactions, wear or particulate debris (such as loose cement, metallic, polyethylene or ceramic particles).
8. Residual pain in the hip.

Intraoperative

1. Femoral shaft or acetabular perforation, fissure, or fracture, which may require the use of internal fixation.
2. Trochanteric fracture.
3. Damage to blood vessels (e.g. iliac, obturator and femoral artery).
4. Temporary or permanent nerve damage (e.g. femoral, obturator or isolated peroneal nerve).
5. Subluxation or dislocation of the hip joint due to implant size or configuration selection, positioning of components and/or muscle and fibrous tissue laxity.
6. Lengthening or shortening of the affected extremity.
7. Surgical delay.

Early Postoperative

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarction.
2. Hematoma and/or delayed wound healing.
3. Subluxation or dislocation.
4. Disassociation of modular components.

Late Postoperative

1. Trochanteric avulsion from excessive muscular tension, weight-bearing, or inadvertent intraoperative weakening of the trochanter.
2. Aggravation of problems in the ipsilateral or contralateral knee and ankle joints due to leg length discrepancy, femoral medialization and/or muscular deficiencies.
3. Acetabular or femoral fracture due to trauma or excessive loading, particularly in the presence of poor bone stock caused by severe osteoporosis, bone defects from previous surgery, intraoperative reaming procedures, or bone resorption.

4. Bone resorption which may contribute to the deterioration of fixation and eventual loosening of the implant.
5. Decrease in joint mobility and range of motion (e.g. maybe caused by periarticular calcification or ossification).
6. Traumatic arthrosis of the ipsilateral knee secondary to intraoperative positioning, of the extremity during surgery.
7. Subluxation or dislocation.
8. Disassociation of modular components.

MRI SAFETY INFORMATION

 MRI Safety Information A patient with the DePuy MARATHON™ XLPE Cemented Cup Hip Prosthesis construct may be scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.	
Name/Identification of device	MARATHON XLPE Cemented Cup Hip Prosthesis
Nominal value(s) of Static Magnetic Field [T]	1.5 T or 3.0 T
Maximum Spatial Field Gradient [T/m and gauss/cm]	30 T/m (3000 gauss/cm)*
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	Whole body transmit coil, Head RF transmit-receive coil
Operating Mode	Normal Operating Mode
Maximum Whole Body	2.0 W/kg (Normal Operating Mode)
Maximum Head SAR	Use of a head coil is not required for scanning of the hip
Limits on Scan Duration	2.0 W/kg whole body average SAR. Scan for up to 15 minutes followed by 15 minutes cooling period (may be repeated up to three times in a single imaging session)

MR Image Artifact	The presence of this implant is not expected to produce an image artifact extending further than 108 mm*
Additional instructions or information essential for safe use in the MR environment.	* Check all components of the hip construct for additional MR labelling information.

DISPOSAL

Safely dispose of all devices used in accordance with local procedures and guidelines.

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